VANNES DE SECTIONNEMENT (GAZ MEDICAUX ET VIDE) SHUT-OFF VALVES (MEDICAL GASES AND VACUUM)

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



VANNE UFM POUR GAZ COMPRIMÉS
UFM VALVE FOR COMPRESSED GASES


VANNE DE COUPURE DE VIDE
VACUUM SHUT-OFF VALVE

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
The shut-off valves are mainly designed to isolate the downstream part of a low-pressure medical gas or vacuum pipeline system.

Other functionalities are usually added to show, locally or remotely, some information about the state of the pipeline system and/or of the valve

6 types of shut-off valves are available and their functionalities are described below.

All the devices described below are medical CE marked.

UFM VALVES:

The UFM assembly is composed of a manual shut-off (which cuts the supply from the source) and of a collector fitted with a pressure gauge (or vacuum gauge in case of vacuum), a pressure transducer and a specific terminal unit which allows the supply of the downstream pipe system in case of problems at the sources level.

VACUUM SHUT-OFF VALVES:

This valve model is fitted with a transducing gauge indicating the actual vacuum level.

AREA SHUT-OFF VALVES:

These valves are a normative and obligatory component of a pipeline system according to the ISO 7396-1 standard.

They enable a physical separation between upstream and downstream without having to dismantle the existing pipes. They are equipped with a specific medical terminal unit which allows the compressed gas supply of the downstream pipe system in case of problems in the upstream portion.

"FIREFIGHTER" VALVES:

These valves usually meet a local regulation of the firefighting services.

It is a manual shut-off valve which enables to know its position remotely (open or closed), and thus the firefighters can operate knowing the state of the gas pipeline system.

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SHUT-OFF VALVES (MEDICAL GASES AND VACUUM)

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BRAS PLAFONNIER / CEILING PENDANT

VALVES FOR RISER AND PRIMARY CIRCUIT:

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They enable to cut the supply of a riser or a primary circuit in order to operate downstream on any problem of the pipeline system.

In view of their location in a pipeline system, they offer a large passage.

CEILING PENDANT VALVES:

The ceiling pendant valves are specific for this application and, in view of their location in a pipeline system, they offer a small passage.

BALL VALVES:

Refer to the specific data sheet.

All these shut-off valves are usually installed in a distribution box of the health facility.

These distribution boxes are sealable, available in several sizes and some are equipped with a screen.

They are prepared and equipped on request.

FEATURES

- Ball valve (Ø from 3/8" to 2 1/2")
- Operating pressure: 10 bar maxi (20 bar maxi from 3/8" to 3/4")
- Collector
- Medical terminal unit (according to model)
- Gauge/Vacuum gauge (according to model)
- Inlet and outlet according to size and type of valve

MATERIAL

- Body and collector: brass and chrome plated brass
- Seals and sealing: according to gas

OPTIONS

- Transducer

MAINTENANCE

Control and maintenance recommendations are detailed in the technical data sheet delivered with the equipment.

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ORDER: Please mention

- Reference
- Gas type or gas mixture components
- Size of the valve
- Selected options and accessories

Example: FIREFIGHTER valve - O₂ - G 1/2" - Transducer

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All views, descriptions and features are given only for information. CAHOUET reserves the right to modify its products without notice.

Distributeur / Distributor

52, rue de Lagny
93100 MONTREUIL - FRANCE
Tél : 33 (0)1.41.72.90.00
Fax : 33 (0)1.41.72.90.01
E-mail : cahouet@cahouet.com
www.cahouet.com